



SCAD-DIC

2D & 3D - DIGITAL IMAGE CORRELATION SYSTEM

FOR NON-CONTACT FULL FIELD DISPLACEMENT,


STRAIN MEASUREMENT & FREQUENCY ANALYSIS.

Abhishek G
Director & CTO
PYRODYNAMICS

Plot No 272 ; Road No 2
KIADB Industrial Area Phase 2 Harohalli;
Taluka Kanakapura; Dist. Ramanagaram
Bengaluru Rural - 562 112


Tel:- +91 9686478833
Email:- pddic.ag@gmail.com
Web:- www.pyrodynamics-india.com
Facebook:- Pyrodynamics
Instagram:- pyrodynamics_india

VIC-3D Calibration Target Grids – Make in India




VIC-3D TEST SETUP

PYRODYNAMICS NOW HAS ESTABLISHED A FACILITY TO MANUFACTURE THE CALIBRATION TARGET GRIDS IN HOUSE THANKS TO MAKE IN INDIA SCHEME INITIATIVE BY GOVERNMENT OF INDIA



VIC-3D - CALIBRATION SCORE COMPARISON :- PD TARGET VS CSI TARGET


Stereo System Calibration



PD - MII - TARGET

Stereo System Calibration

Target 10mm (H95-00) + Edit Automatically recalibrate Recalibrate



Imaging system Standard, radial 2 Coordinate system Auto

Analyze Report Accept Cancel

CSI - TARGET

1	2	^
51	0.031	0.032
52	0.045	0.043
53	0.048	0.045
54	0.032	0.036
55	0.015	0.023
56	0.030	0.035

Value Result Std ^

Camera 0

Type	Pinhole
Center (X)	1196.26
Center (Y)	1213.30
Focal Le...	5127.99
Focal Le...	5127.76
Avg. ma...	8.09692
Min. ma...	7.91201
Max. m...	8.48167

Camera 1

< >

Score: 0.026

1	2	^
36	0.017	0.028
37	0.025	0.021
38	0.025	0.023
39	0.017	0.023
40	0.026	0.035
41	0.026	0.027

Value Result Std ^

Camera 0


Type	Pinhole
Center (X)	1157.35
Center (Y)	1261.60
Focal Le...	5140.32
Focal Le...	5142.73
Avg. ma...	8.03121
Min. ma...	7.82789
Max. m...	8.39439


Camera 1

< >

Score: 0.029

VIC-3D RESULTS - PD CAL TARGET V/S CSI CAL TARGET





**Designed & Manufactured
by Pyrodynamics**



Low Speed Camera Mounting



High Speed Camera Mounting

(Image Courtesy:- Kuwait University – Kuwait)


DIC Measures Displacements.


Strain is a Derivative of Displacement. Lagrange Strain tensor is used for calculation of Strains

Calculation of strains:


$$\epsilon_x = \frac{\partial u}{\partial x}, \quad \epsilon_y = \frac{\partial v}{\partial y},$$

$$\gamma_{xy} = \frac{\partial u}{\partial y} + \frac{\partial v}{\partial x}$$





Two USB CMOS Cameras
Interfaced to the USB Port of the
PC.




Displacement and Strain
Measurement In-Plane
Direction
(X and Y Direction)
And
Out of Plane Direction
(Z Direction)

Results Obtained from 3D Digital Image Correlation

Contour Co Ordinates:- X, Y and Z ; Displacement:- u, v and w.

Strains:- ε_{xx} , ε_{yy} , ε_{xy} , ε_1 , ε_2 , Von Mises, Tresca

and Directions of Principal Strains



Fire Camera interfaced to the
Fire Wire Port of the PC.

Displacement and Strain Measurement
In-Plane Direction
(X and Y Direction).

Results Obtained from 2D Digital Image Correlation

Contour Co Ordinates :- X and Y

Displacement:- u and v.

Strains:- ϵ_{xx} , ϵ_{yy} , ϵ_{xy} , ϵ_1 , ϵ_2 Von Mises, Tresca and

Directions of Principal Strains

List of Available CMOS USB 3.0 Cameras



Specification	CSI-04-0200	CSI-04-0202	CSI-04-0203	CSI-04-0204
Camera Interface	USB3.0	USB3.0	USB3.0	USB3.0
Camera Chip	CMOS	CMOS	CMOS	CMOS
Camera Sensor	Sony IMX174	Sony IMX250	Sony IMX255	Sony IMX253
Maximum Resolution	2.3 MP - 1920 X 1200	5.0 MP - 2448 X 2048	8.9 MP - 4096 X 2168	12.3 MP - 4096 X 3000
Maximum Frame Rate	164 fps	75 fps	43 fps	30 fps
Pixel size	5.86 µm	3.45 µm	3.45 µm	3.45 µm
Sensor size	1/1.2"	2/3"	1"	1.1"
Lens Mount	C	C	C	C
Compatible with Active Optic Cables	Yes	Yes	Yes	Yes



Correlated Solutions Inc GigE Cameras for VIC-3D, 3D DIC System.

These Cameras operate only with Desktop PC's and a 2 Port GigE PCIe Interface Card needs to be installed in the PC.

Specification	CSI-04-0235	CSI-04-0236	CSI-04-0237	CSI-04-0238
Camera Interface	GigE	GigE	GigE	GigE
Camera Chip	CMOS	CMOS	CMOS	CMOS
Camera Sensor	Sony IMX250	Sony IMX255	Sony IMX254	Sony IMX342
Maximum Resolution	2448 x 2048 (5 MP)	4096 X 2160 (8.9 MP)	4096 X 3000 (12.3 MP)	6464 X 4852 (31 MP)
Maximum Frame Rate	162 fps	93 fps	68 fps	26 fps
Pixel Size	3.45 µm	3.45 µm	3.45 µm	3.45 µm
Sensor Size	2/3"	1"	1.1"	APS-C
Lens Mount	C	C	C	F




Correlated Solutions Inc High Speed Cameras for 2D/3D DIC Systems.

Specification	CSI-04-0188 Photron AX50 170K-M-16GB	CSI-04-0178 Photron AX100 540K-M-16GB	CSI-04-0061 Photron UX50 160K-M-16GB	CSI-04-0064 Photron UX100 800K-M-16GB
Camera Interface	Ethernet	Ethernet	Ethernet	Ethernet
Camera Chip	CMOS	CMOS	CMOS	CMOS
Maximum Resolution	1024 x 1024	1024 x 1024	1280 x 1024	1280 x 1024
Maximum Frame Rate	2000 fps	4000 fps	2000 fps	4000 fps
Lens Mount	F	F	F	F
Compatible with Active Optic Cables	No	No	No	No

Technical Data



LED Video Light	SL150II	SL200II
Power Supply	AC100V-240V (50/60Hz)	
Max Power	150W	200W
Color Temperature	5600±200K	
100% Light Brightness(LUX,1M)	≈58000	≈74000
CRI ≈	96	
TLCI ≈	97	
Light Brightness Range	10%-100%	
Operation Temperature	-10°C ~ 40°C	
Channel	32	
Group	16(A, B, C, D, E, F, 0-9)	
Dimension (not included a standard reflector)	34x20x16cm	36.7x20x16cm
Net Weight	2.95Kg	3.33Kg



Remote Control	RC-A6
Power Supply	3.0V(AAA*2)(not include)
Wireless Frequency	2.4GHz
Distance (emptiness)	50m
Channel	32
Group	16 groups(A-F, 0-9)
Standby Time	Over one year
Operation Temperature	-10°C ~ 40°C
Dimension	120*38*15mm
Net Weight	30g



50mm Lens



85mm Lens



100mm Lens



C to F Mount Adapter



135mm Lens



Lets Think Collectively

Thank You

Abhishek G
Director & CTO
PYRODYNAMICS
Plot No 272
KIADB Industrial Area Phase 2
Harohalli ; Taluka Kanakapura
Dist. Ramanagaram
Bengaluru Rural - 562 112
Karnataka State – India

Tel:- +91 9686478833

Email:- pddic.ag@gmail.com
LinkedIn & Facebook:- Pyrodynamics

**Moto of
Pyrodynamics**

**“A Company of
Value Rather
than a Company
of Success”**